



The Twinkle Factory

Stain different, tag FAST.

greenFAST

GreenFAST is a small protein tag of 125 amino acids (14 kDa) engineered from the photoactive yellow protein (PYP) from *Halorhodospira halophila*. GreenFAST can be fused to the N- or C-terminus of any protein of interest. greenFAST forms non-covalent and reversible fluorescent molecular assemblies with ^{Tf}Lime, enabling to fluorescently label greenFAST-tagged proteins in vitro, in live cells and in multicellular organisms.

greenFAST sequence

MEHVAFGSEDIENLAKMDDEQLDGLAFGAIQLDGDGNILQYNAAEGDITGRDPKQVIGKNFFKDVATGTDSPFYRKFKKE
GVASGNLNTMFEWMIPTSRGPTKVKVHMKKALSGDSYWVFKRV

DNA sequence coding for greenFAST

atggagcatgttgcccttggcagtgaggacatcgagaacactctggccaaaatggacgacgaacaactggatgggttggccttggcgcaattcagctcga
tggcgacgggaatatcctgcagtacaatgctgctgaaggagacatcacaggcagagatcccaaacaggtgattgggaagaacttctcaaggatgttga
actggaacggattctcccagttttaccgcaaattcaaggaaggcgtagcgtcagggaaatctgaacacatgttcgaatggatgataccgacaagcaggg
gaccaaccaaggtcaaggtgcacatgaagaaagcccttccggtgacagctattgggtcttgtgaaacgggtg

Notice to Buyer/User

The Buyer/User has a non-exclusive license to use this system or any component thereof for research use only. Commercial use of this system or any components thereof requires a license from Twinkle Bioscience S.A.S, 45 rue d'Ulm, 75005 Paris, France. For detailed information, e-mail contact@the-twinkle-factory.com. The products and/or their use may be covered by one or more of the following patents and patent applications:

- EP 3 164 411; JP 2017-527,261; WO 2016/001,437; US 10,138,278 (Fluorogen activating and shifting tag (FAST))
- EP 3 404 022; JP 2020-520,676; WO 2018/211,090; US 2020-124,611 (Membrane-impermeant fluorogenic chromophores)
- EP 3 719 007 (Split photoactive yellow protein complementation system and uses thereof)

These patents and patent applications are owned by CNRS (France), PSL (France), Sorbonne Université (France), and/or Institut Curie (France), and exclusively licensed to Twinkle Bioscience.